

# STERISHEET STERILIZATION WRAPS

## STERISHEET 340 SMS



### DESCRIPTION

Sterisheet sterilization wraps are best in class Sterile Barrier Systems for CSSDs in hospitals and clinics. Sterisheet SMS products are available as interleaved, bonded, or single sheets.


### COMPOSITION

100% synthetic fibers made of polypropylene.

Spunbond / Meltblown / Spunbond

Spunbond network consists of long, strong, and thick polypropylene filaments.

Meltblown layers are made of short and thin polypropylene microfibers.



### SUITABLE FOR STERILIZATION METHODS

- Steam
- EO
- Vaporized hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)

### APPLICATION

For small and medium-size trays, light packs, instrument kits.

### SIZES AVAILABLE

Choose standard sizes to optimize your costs.

Other sizes are also available upon request.

75x75 cm            90x90 cm

100x100 cm        120x120 cm

# STERISHEET STERILIZATION WRAPS

## STERISHEET 340 SMS



### COLORS



### PERFORMANCES

Mechanical and bacterial barrier properties are internally tested on a routine basis according to regulatory requirements of all of our products. Combining them with random tests conducted by external and accredited laboratories leads Sterimed to secure the best performances on Sterisheet products.

#### Mechanical Properties

Preservation of pack integrity from closure till the point of use depends on the material's resistance to tearing, puncturing, and breaching stresses generated all along the distribution with the hospital. Any mechanical weaknesses will increase the risks of event-related ingress of microorganisms into the pack. Excellent mechanical properties will provide you additional safety while using our materials. Optimal strength and resistance provided in every sheet.

PROPERTIES	TYPICAL
Basis Weight	40 g/m <sup>2</sup>
Air permeability	150 l/min/dm <sup>2</sup>
Thickness	300 µm
Hydrost Test	40 mBar
Tensile strength MD	2.0 kN/m
Tensile strength CD	0.9 kN/m
Elongation MD	65 %
Elongation CD	65 %
Burst	200 kPa
Tears SM	3200 mN
Tears ST	5800 mN

#### Bacterial Properties

Sterilization wraps must prevent microorganisms' ingress inside the package. To reach this performance, different types of testing have to be performed on the products to reproduce both kinds of ingress vehicles possible:

- Airborne ingress
- Waterborne ingress

For instance, according to the TNO final pack test method, validating the materials & the folding technique, STERISHEET 340 SMS exhibits a BARRIER PROTECTION > 99,99%

# STERISHEET STERILIZATION WRAPS

## STERISHEET 340 SMS



### COMPLIANCE TO STANDARDS

Sterisheet products range is classified as a Class I Medical Device according to the European Medical Regulation (MDR 2017/745/EU). Its CE marking illustrates the relevant compliance.

Sterisheet SMS products conform with the standards below:

EN ISO 11607-1

EN 868-2

### OUR WRAPS MANUFACTURING CERTIFICATION

ISO 13485 standard

### PACKAGING PRIOR TO USE

Sheets displays is optimized by adjusted folding depending on the size and type of the product. We have carefully tested the best solution to ensure the most convenient handling for end users.

- PRIMARY TRANSPORT PACKAGING

Number of sheets is optimized and wrapped in a transparent polyethylene bag with quick product ID

- SECONDARY TRANSPORT PACKAGING

Secondary packaging is a neutral brown color cardboard box providing transportation stress resistance

### LABELING

Product traceability is fully insured through labeling according to regulations on each transport packaging.

### STORAGE CONDITIONS

Sterimed recommends the following storage conditions for best performance of sterilization wraps: Storage in a cool, dry location away from direct exposure to natural light, strong artificial light & UV sources. Cardboard boxes should never be stored in direct contact with the floor. Storage of the products shall be done in areas that are not subject to extreme temperature changes such as in contact with heated objects, vents, or cold walls.

# STERISHEET STERILIZATION WRAPS

## STERISHEET 340 SMS



As per AAMI ST79 “*Comprehensive guide to steam sterilization and sterility assurance in health care facilities*” recommendations: Before use, hold packaging materials at room temperature (20°C to 23°C) and at a relative humidity ranging from 30% to 60% for a minimum of 2 hours is a good practice for optimum use performances.

### USE BY DATE

Provided the above storage conditions are met, the upper limit of the time interval during which the performance characteristics of the sterilization wrap are demonstrated is 5 years of the manufacturing date.

### ENVIRONMENTAL IMPACT & WASTE MANAGEMENT

Oil derived product, renewable content 0%.

Disposal as per local regulations after use.